

Cat. No.	Descrizione	Pannello	no. Test	Pack.
213043	Elettroliti	CO ₂ , K ⁺ , Na ⁺ , CA ₂ ⁺ , P, Mg ₂ ⁺ , Cl ⁻	7	10
213077	Chimica generale 4	TP, ALB, GLO, A/G, TBIL, DBIL, IBIL, ALT, AST, GGT, ALP, UREA, CRE, UA, TG, CHOL, GLU, HDL-C, LDL-C	19	10
213082	Glucosio e Lipidi	TG, CHOL, HDL-C, LDL-C, GLU, GSP	6	10
213078	Funzione epatica	TP, ALB, GLO, A/G, TBIL, DBIL, IBIL, ALT, AST, GGT, ALP	11	10
213079	Funzione renale - su richiesta -	ALB, CRE, UREA, UA, CO ₂ , Ca ₂ ⁺ , P	7	10

TP	Proteine Totali
ALB	albumina
GLO	globulina
A/G	rapporto albumina/globuline
TBIL	bilirubina totale
DBIL	bilirubina diretta
IBIL	bilirubina indiretta
ALT	transaminasi
AST	aspartato transaminasi
GGT	gamma GT
ALP	fosfatasi alcalina
TG	trigliceridi
CHOL	colesterolo
HDL-C	lipoproteine alta densità
LDL-C	lipoproteine bassa densità
GLU	glucosio
GSP	Fruttosamina
CRE	creatinina
UREA	azotemia
UA	Acido urico
CO ₂	anidride carbonica
CA ₂ ⁺	Calcio
P	Fosforo
K ⁺	Potassio
Na ⁺	Sodio
Mg ₂ ⁺	Magnesio
Cl ⁻	Cloro